


COVID-19 Self-Testing Instructions BTNX Rapid Response™

Follow these six steps to complete the COVID-19 self-test using the BTNX Rapid Response™ rapid antigen test kit.

Step 1. Prepare supplies



Wash hands before starting. Prepare your testing area and lay out the supplies you need on a clean surface: 1. Test tube and cover 2. Buffer fluid in ampule 3. Test device 4. Tube rack 5. Nasal swab

Step 2: Prepare test tube




Empty one ampule of buffer fluid into the test tube. Place tube in holder.

Step 3. Obtain Throat and Nares Swab



Swab the back of the throat (avoiding the tongue, teeth and lips) in the tonsillar area for 5 seconds on each side. Swab the nares using the same swab. Insert swab about 1 cm into the nares, rotate 5 times, remove and repeat on the other side.

Step 4. Place the swab into the tube




Place the swab into the tube and stir the swab around in the fluid. Squeeze the swab through the tube 10-15 times. Place in the holder for 2 minutes. Remove nasal swab from tube while squeezing the fluid out. Place the cap firmly on the tube.

Step 5. Apply fluid to test device



Squeeze **3 drops** of fluid from the test tube onto the sample well of the test device.

Step 6. Read test result



Read results after 15 minutes. One red line (C) - indicates a negative result, Two red lines (C and T) - indicate a positive result. One line (T) - is an invalid test (re-do the test).

Positive test? Positive results do not require a confirmatory PCR test. Refer to the [Ministry of Health's Self-Testing Guidelines](#) for advice on how to proceed once results are displayed and test processing time is complete.